




RelaxSan Cotton Socks for Man is a line of cotton socks for man with a graduated compression, for the prevention and the treatment of problems relating to legs blood circulation.

RelaxSan items have been originated after years of experience and are produced following the recommendations of the European community rules for the Class I. Any items are Certified CE.

RelaxSan items are produced with yarns of very high quality, to guarantee a long duration and have an anallergic action at skin contact.

Sanitized treatment guarantees an antibacterial action constant in time.



Available Items:

RelaxSan Cotton Socks 18-22 mmHg

RelaxSan Cotton Socks 22-27 mmHg

RelaxSan Cotton Socks 18-22 mmHg X-Static



The item with X-Static Silver yarn, besides graduated compression, has also the properties of silver, listed here below:



NATURAL

100%. It's safe, non toxic, does not contain any chemical or pesticides residues.



ODOUR CONTROL

Neutralizes ammonia, denatured proteins and blocks bacterial proliferation.



THERMO-REGULATOR

Fresh in summer and warm in winter for its thermal conductivity; both the reflecting and the radiating power improve humidity passage.



THERAPEUTIC

Thanks to its conductive characteristics this new fibre gives some benefits to health and improves blood circulation.



ANTIMICROBIAL

Removes those bacteria which causes smell and removes fungi responsible for the "athlete's foot".



ANTISTATIC

For its high electric conductivity.

RelaxSan Cotton Socks items, thanks to their graduated compression, improve blood circulation and are recommended for those people who make sedentary activities and travel frequently.



Compression categories:



18-22 mmHg* Compression recommended for the therapy of the light chronic venous insufficiency, varicose veins and swelling.



22-27 mmHg* Compression recommended for the therapy of the venous insufficiency, varicose veins, oedema. Phlebitis, after surgeries and in the prevention of the Thrombosis.

* The compression recommended refers to ankle, and then decreases gradually.

How to wear the Cotton Socks:



- 1 Delicately, enrol the item from the elastic to the foot tip.
- 2 Position correctly the tip and the heel of the socks on the foot.
- 3 Accurately and delicately, spread out the item from the ankle to the knee.
- 4 Distribute evenly the item removing any crease.

Sizes Table

Tg - Size	2°	3°	4°	5°	6°
European	37-39	39 - 41	41 - 43	43 - 45	45 - 47
American	5 - 7	7 - 8	8 - 10	10 - 12	12 - 14
British	4 - 6	6 - 7	7 - 9	9 - 11	11 - 13
Asian	23,5-25,5	25,5-26,5	26,5-28	28 - 30	30 - 32

SIZE	ANKLE CIRCUMFERENCE
2	18-21
3	21-24
4	24-27
5	27-30
6	30-33



by

GT CALZE s.r.l. - Via W. Tobagi 17/21 - 46040 Casaloldo (MN) MADE IN ITALY

